



QualityPatrol

Quality Plans

Specifications

Configure

QualityPatrol

Configure

Archiving

Define the configuration settings for the PatrolSuite Archiving Engine (separate license required).

Operations

Define the steps of the manufacturing process where testing is performed.

Quality States

Define the rules for how samples can be processed.

Quality Plans & Specs

Define what tests need to be completed when a raw material is received or a product is made, and the corresponding specifications.

[Quality Plans & Specs](#)
[Quality Plan Statuses](#)

QualityPatrol

Quality Plans

Search quality plans

Search

Reset

StatusActive

FacilityAny

Product GroupAny

OperationAny

ProductEnter name or number...

Test MethodEnter name or number...

Business PartnerEnter name or number...

123456...874

Page

Go

+ New Quality Plan

Product	Operation	Business Partner
12346 Better Butter Crème Deluxe TCCS-BETTERCREME	ARL-BIB	

Quality Plan

12346 ARL-BIB Active

HistoryEditCloneDelete

Product12346 Better Butter Crème Deluxe

OperationARL-BIB

CustomerAny

Shelf Life

Quality Plan StatusApproved - Commercialized

Comments
test comment

Related Quality Plans0

Test Methods

Specifications

Frequencies

Locations

Monitoring

Sources

Calculations

Rules

Quality Plan Calculations

Σ VS 11113

=

VM 0003 Average



Contents

CHAPTER 1: INTRODUCTION..... 4

NEW INTRODUCTION..... 4

 ACCESS 4

 BROWSER 4

NEW DEFINITIONS..... 4

 SYMBOL KEY 8

CHAPTER 2: SECURITY 9

 INTRODUCTION..... 9

CHAPTER 3: QUALITY PLAN SEARCH 10

NEW INTRODUCTION..... 10

 VIEW 10

 SEARCH 11

CHAPTER 4: QUALITY PLAN OVERVIEW 12

NEW QUALITY PLANS - OVERVIEW 12

CHAPTER 5: QUALITY PLAN - CREATE NEW..... 13

NEW CREATE NEW QUALITY PLAN 13

CHAPTER 6: QUALITY PLAN - EDIT 15

NEW HISTORY..... 15

NEW EDIT QUALITY PLAN INFORMATION 15

NEW CLONE QUALITY PLAN 16

NEW DELETE QUALITY PLAN 18

CHAPTER 7: RELATED QUALITY PLANS..... 19

NEW VIEW 19

CHAPTER 8: TEST METHODS..... 20

NEW ADD TEST METHODS..... 20

CHAPTER 9: SPECIFICATIONS..... 23

NEW SPECIFICATION LEVELS 23

NEW SPECIFICATION TYPES 23

NEW SPECIFICATION TOLERANCE LIMITS 24

NEW TEST METHOD SPECIFICATIONS 24

NEW SPECIFICATIONS AUDIT TRAIL 27

CHAPTER 10: SPECIFICATION EDITOR 29

NEW INTRODUCTION..... 29

NEW VIEW 29

NEW DEFINE NEW TEST METHOD SPEC..... 30

NEW EDIT EXISTING TEST METHOD SPEC..... 32

NEW DEFINE NEW PRODUCT SPEC..... 33

QualityPatrol – Quality Plans & Specs

User Guide

Version 2025-10



NEW	EDIT EXISTING PRODUCT SPEC	34
NEW	DEFINE NEW CUSTOMER SPEC	36
NEW	EDIT EXISTING CUSTOMER SPEC	38
CHAPTER 11: SEARCH SPECIFICATIONS		40
NEW	INTRODUCTION	40
NEW	SEARCH	40
NEW	EDIT	41
CHAPTER 12: SAMPLE FREQUENCIES		42
NEW	INTRODUCTION	42
NEW	VIEW	42
NEW	STARTUP	42
NEW	FREQUENCY SELECTION	43
CHAPTER 13: TEST LOCATIONS		45
NEW	INTRODUCTION	45
NEW	VIEW	45
NEW	EDIT	45
CHAPTER 14: MONITORING		47
NEW	INTRODUCTION	47
NEW	VIEW	47
NEW	FACILITY OVERVIEW	47
NEW	LINE DETAILS	48
CHAPTER 15: SOURCES		49
NEW	INTRODUCTION	49
NEW	VIEW & EDIT	49
CHAPTER 16: CALCULATIONS		51
NEW	INTRODUCTION	51
NEW	CREATE EXPRESSION	51
NEW	TEST EXPRESSION	53
CHAPTER 17: QUALITY STATE RULES		55
NEW	INTRODUCTION	55
NEW	VIEW	56
NEW	NEW RULE	56



Chapter 1: Introduction

NEW Introduction


QualityPatrol™ is a PatrolSuite™ module that facilitates:

- The definition of the product quality plan, including process conditions and product testing requirements and the corresponding targets and specifications. In regulated industries such as medical devices, this is the design transfer step from the device master record (DMR) to manufacturing.
- The capture and storage of the information required by the quality plan (e.g. the test results) and the subsequent determination of the product's quality state. In regulated industries such as medical devices, this is the device history record (DHR).

This user guide covers the creation of a product **quality plan** and the corresponding **specifications**. The remainder of the QualityPatrol configuration settings is covered in the “QualityPatrol – Configure” user guide.

Access

To access the PatrolSuite homepage, obtain the URL and login credentials from your IT department.

 PatrolSuite module access is governed by both PatrolSuite and Active Directory security protocols. Therefore, users will see different homepage content and will have different access levels based on assigned roles. For assistance with security or permissions, contact your IT department.

Browser

PatrolSuite™ is accessed through a web browser. The platform and all its applications are certified for use with Google Chrome, Microsoft Edge, and Safari. To ensure full functionality and the best user experience, MAP strongly recommends using one of these supported browsers.

NEW Definitions

Business Partner – Customers and Suppliers.

Business Unit – Grouping of manufacturing facilities according to the market to which their products are sold.

QualityPatrol – Quality Plans & Specs

User Guide

Version 2025-10



Comprehensive Testing – If an administrator has turned comprehensive testing “on” for a given test location, then all test results must be recorded before the user is allowed to save/close the edit test results window. This feature can be used to help ensure that all tests are performed as required.

Cost of Quality (CoQ) – The total cost of ensuring and maintaining product quality. It includes all costs involved in preventing defects, appraising product quality, and dealing with failures. CoQ is divided into four categories:

1. Prevention Costs (Quality Assurance)
 - Costs to prevent defects before they happen such as training, process design, quality planning, policy and procedure documentation, and preventive maintenance.
2. Appraisal Costs (Quality Control)
 - Costs of measuring and monitoring product quality such as product inspection and testing, audits, calibration of instruments, statistical process control (SPC).
3. Internal Failure Costs (Nonconforming Product)
 - Costs of defects **found before** the product reaches the customer such as scrap, rework, downtime, re-inspection.
4. External Failure Costs (Customer Complaints)
 - Costs of defects **found after** the product is delivered to the customer such as customer complaints, warranty claims, returns, product recalls, loss of reputation.

Cost of Poor Quality (CoPQ) – Avoidable costs that result from delivering a substandard product or service. CoPQ = Internal Failure Costs + External Failure Costs. Accordingly, CoPQ is a subset of the CoQ often referred to as “the cost incurred when things go wrong.”

Currency – A system of money in common use within a particular country or economic region, used as a medium of exchange, store of value, and unit of account in financial transactions.

Defect – Flaws, faults, or deviations in a manufactured item that prevent it from meeting design specifications, quality standards, or customer expectations.

Defect Category – Grouping of defects to help identify, analyze, and correct issues systematically. Common defect categories include:

- Aesthetic or Cosmetic Defects - Visual imperfections that do not affect the product’s functionality but can impact customer perception and marketability.
- Design Defects - Flaws inherent in the product's design that make it unsafe, ineffective, or prone to failure, even if manufactured perfectly.
- Documentation/Labeling Defects - Inaccurate or missing product information, which may lead to misuse, safety issues, or regulatory non-compliance.
- Functional Defects - The product does not perform as intended or fails under specific conditions.
- Manufacturing Defects - Flaws introduced during the production or assembly process, even when the design is correct.
- Material Defects - Deficiencies in raw materials or components that affect product performance, safety, or durability.
- Packaging Defects - Issues related to how the product is packaged for shipping, storage, or display.
- Process Defects - Defects caused by errors in the manufacturing process settings, procedures, or conditions.

Department – An organized unit within a company that is responsible for managing a specific set of functions, processes, or activities that contribute to the organization’s overall goals and operations.

QualityPatrol – Quality Plans & Specs

User Guide

Version 2025-10



Disposition – The remediation instructions assigned to a product with a rejected quality state. The goal of the remediation is to minimize the financial loss due to making product that does not meet quality requirements.

Enterprise Resource Planning (ERP) – Software system that enables management of the day-to-day business activities such as accounting/finance, procurement, supply chain, human resources, sales order conversion to manufacturing orders, distribution, project management, compliance, and others.

Facility – A physical location where raw materials are converted into finished or semi-finished products through various production processes, machinery, labor, and technology.

Manufacturing Execution System (MES) – Software system that tracks and documents the transformation of raw materials into finished goods, including inventory consumption and product genealogy. In a world class manufacturing environment, the MES is integrated to communicate bi-directionally with the ERP and the Quality System (PatrolSuite).

Market (“Customer Market”) – Groups of consumers that a company targets to sell its products or services, based on shared needs, characteristics, or behaviors. These markets help companies focus their marketing, product development, and sales strategies by understanding who their customers are and what they need.

Operation – Location where a sample originates for which testing is performed.

Order – A set of data with a unique order number (typically originating from the MES or ERP) to identify a:

- Inbound Receiving order (often a “purchase order”) for raw materials
- Manufacturing order and its related product and customer
- Outbound Shipping order for finished goods being shipped to customers

Product – A tangible output or item that is created through a controlled process involving raw materials, labor, machinery, and other inputs, intended for use, sale, or further processing. Types of products include:

- Finished Goods: Ready for sale to end users (e.g. computers, bicycles, paper).
- Semi-Finished Goods: Used as components in other products (e.g. engine parts, circuit boards).
- Raw Materials: Basic inputs for production (e.g., glue, steel, plastic pellets, lumber).

Product Group – A collection of related products that share common characteristics, functions, markets, manufacturing processes, or branding, and are managed together for strategic, operational, or marketing purposes.

Quality Assurance (QA) – A proactive process that focuses on preventing defects by ensuring that the processes used to manage and create deliverables are effective and followed correctly.

- Focus: Process-oriented
- Goal: Prevent defects before they happen
- When: Throughout the product development and prior to the manufacturing process
- Methods: Process audits, training, standard operating procedures (SOPs), continuous improvement
- Responsibility: Business management, quality management, and process engineers

Quality Control (QC) – A reactive process that focuses on identifying defects in the finished product through inspection and testing.

- Focus: Product-oriented
- Goal: Detect and fix defects after they occur
- When: After production or at specific checkpoints during production
- Methods: Inspections, measurements, product testing, sampling, statistical process control (SPC)
- Responsibility: Laboratory technicians, quality inspectors, or testing personnel

QualityPatrol – Quality Plans & Specs

User Guide

Version 2025-10



Quality Plan – A formal set of electronic files or settings that define the specific quality practices, resources, standards, procedures, responsibilities, and inspection or testing methods to be applied to a particular product or process during manufacturing.

Quality State – The assessment of a sample's conformance to quality requirements. For example, a sample for which all tests are completed, and all are in spec, might have a quality state of "Approved."

Raw Material – A basic, unprocessed, or minimally processed substance that is used as the starting input in the production or manufacturing of goods and products.

Reasonable Data Limit (RDL) – The range in which all test results would normally be expected to fall. If an RDL is defined, and a test result is entered that exceeds the RDL, the user will get a warning message.

Sample – A QualityPatrol record reflecting a point in time at which product characteristics (i.e. test results) or process parameters (i.e. process conditions) are recorded. Samples are related to the Operation and Workcenter where they are being manufactured.

Sample Status "Incomplete" – A sample that has one or more data yet to be collected. For example, if a Sample requires 5 tests to be performed, but only 4 tests have been completed, then the sample is incomplete.

Sample Status "Completed" – A sample for which all tests are completed. Note that complete/incomplete has nothing to do with whether the results are in spec or out of spec.

Sample Location – A physical location on a single sample from which tests are performed for X-bar testing.

Sample Type "Start Up" – Identifies a sample of product that is not intended to be sold to a customer or used in further processing. **Use Case:** When a manufacturing line begins producing "Product 123," the first few samples are identified as Start Up and are used to determine when the product meets quality requirements. Once quality requirements are met, then the next sample is identified with the sample type "Production." **Note:** Start Up samples can be and are typically omitted from statistical analysis (DataPatrol & DashPatrol).

Sample Type "Production" – Identifies a sample of product that is expected to meet quality requirements and be sold to a customer or used in further processing.

Specification Limit – The product measurements within which confirm the acceptability of a product, or the process conditions within which likely result in the production of acceptable product. Product specification limits are typically defined by the customer. Process specification limits are typically defined by the manufacturer.

Specification Limit "Upper" – The highest value that a test result can be and still be considered acceptable to the customer.

Specification Limit "Lower" – The lowest value that a test result can be and still be considered acceptable to the customer.

Test Location – A physical location within the Facility where tests are completed.

Unit – The smallest measurable or sellable quantity of a product. Units are typically used for production, inventory, pricing, sales, and distribution purposes.

Unit of Measure – A standard quantitative unit used to specify, track, and manage the amount of a product,

QualityPatrol – Quality Plans & Specs

User Guide

Version 2025-10







material, or resource in the manufacturing, inventory, sales, or procurement processes.

Warning Limit – A numerical value that, when exceeded, raises a flag prompting a timely, proactive response to eliminate or minimize the possibility of the parameter from reaching an alarm limit such as a specification.

Workcenter – A sub-division of an Operation. An Operation can have one or more Workcenters. Often, a workcenter reflects a unique manufacturing machine or asset.

Symbol Key

The following symbols are used in this user guide.

Typeface	Description
	Indicates a “required” data entry field within a UI.
	Indicates a new feature, or new information regarding existing functionality.
	Indicates rules regarding access security.
	Click the help icon to view additional information.



Chapter 2: Security

Introduction

See the “QualityPatrol – Configure” user guide for information regarding QualityPatrol security.



Chapter 3: Quality Plan Search

NEW Introduction

For each **Product** (or raw material) and its associated **Operation**, a **Quality Plan** defines the following:

Required Tests - Specifies which quality tests must be conducted to ensure compliance with product and process standards.

Timing of Tests - Indicates when each test must be performed (e.g., before processing, during production, post-production, at receiving).


Location of Tests - Specifies where each test should take place, such as:

- At the point of manufacture (on-line or in-process)
- In a laboratory
- At receiving or final inspection areas

Applicable Specifications - Lists the standards or specifications (e.g., ASTM, ISO, internal company specs) that the product or material must meet.

View

1. To view the existing quality plans, click on the Quality Plans hyperlink in the QualityPatrol card on the PatrolSuite homepage.



The quality plan search page will appear.

QualityPatrol

[Quality Plans](#)

Search quality plans

Status	Active	Facility	Any
Product Group	Any	Operation	Any
Product	Enter name or number...	Test Method	Enter name or number...
Business Partner	Enter name or number...		



Search

- To find a quality plan of interest, click on the Quality Plans hyperlink in the QualityPatrol card on the PatrolSuite homepage, use the filters, or the Search field and click on the 'Search' button.

QualityPatrol
Quality Plans

Search quality plans [Search] [Reset]

Status: Active Facility: Any

Product Group: Any Operation: Any

Product: Enter name or number... Test Method: Enter name or number...

Business Partner: Enter name or number...

Note: The Product, Business Partner, and Test Method fields are partial search fields.

Note: The Test Method field will return all quality plans that contain the entered test method.

- Click on the quality plan name hyperlink to open the quality plan.

QualityPatrol
Quality Plans

Search quality plans [Search] [Reset]

Status: Active Facility: Any

Product Group: Any Operation: Any

Product: Enter name or number... Test Method: Enter name or number...

Business Partner: Enter name or number...

1 2 3 4 5 6 7 ... 874 Page Go [New Quality Plan]

Product	Operation	Business Partner
12346 - Better Butter Crème Deluxe TCCS-BETTERCREME	ARL-BIB	



Chapter 4: Quality Plan Overview

NEW Quality Plans - Overview

5. A quality plan is organized into four areas:

Breadcrumb Menu. A navigational aid used to show users their current location within a hierarchical structure of a website. Helps users understand the path they have taken and how to navigate back to previous pages.

Quality Plan Information = Associates the quality plan with a product, operation, and, if applicable, a customer.

Related Quality Plans. Shows plans for the same product, but different operations or customers.

Quality Plan Tabs = Defines the applicable test methods, specifications, and corresponding configuration settings.

Quality Plans > 04810 > ARL-BIB

Quality Plan

04810 > ARL-BIB Active

History Edit Clone Delete

Product

04810 ▾ HOC CHP CDO #240/1Z

Quality Plan Status

None

Operation

ARL-BIB

Comments

Customer

Any

Shelf Life

▼ Related Quality Plans 3

ARL-BIB

04810 This Plan

MUR-Cookie-Mixer

04810

MUR-Cookie-Pkg

04810

MUR-Cookie-Table

04810

Test Methods Specifications Frequencies Locations Monitoring Sources Calculations Rules

Test Type All ▾ Product Specifications Add Test Method Edit


Test Method	Certificate	Control Class	Active
VS 2100-Fats-CEM ▾ Fats - CEM SmartTrac	■	Inherit [Standard]	<input checked="" type="checkbox"/>




Chapter 5: Quality Plan - Create New

NEW Create New Quality Plan

6. To create a new quality plan, click on the Quality Plans hyperlink in the QualityPatrol card on the PatrolSuite homepage.



- a. Click on the 'New Quality Plan' button.



- b. Add the required information and click the 'Save Changes' button.

Note: Click into the Product field and select the applicable product from the product picker.

Note: Select the applicable Operation from the dropdown.

Note: Click into the Customer field and select the applicable customer from the customer picker. If the product, as tested via the quality plan and corresponding specifications, can be shipped to any customer, then leave the Customer field blank.

Note: Enter applicable Quality Plan Comments. These comments are visible to the user performing testing on the sample.

Note: The quality plan status is typically first set to "Do not use – Initial Draft" so that it does not get consumed before it has been completed. The status can be changed at any time while managing the details of the plan.

QualityPatrol – Quality Plans & Specs

User Guide

Version 2025-10



Quality Plans > New Quality Plan

Quality Plan

New Quality Plan ☒ Active

Product
703335 Cheesecake NY Strawberry

Operation
MOR-P2-Cake5-Pkg

Customer
0070085393 DJVM LLC

Quality Plan Status
1. Do Not Use - Initial Draft

Comments

The new quality plan will be displayed with the Test Methods tab selected.

Quality Plans > 703335 • MOR-P2-Cake5-Pkg • 0070085393

Quality Plan

703335 • MOR-P2-Cake5-Pkg • 0070085393 Active

Product
703335 Cheesecake NY Strawberry

Operation
MOR-P2-Cake5-Pkg

Customer
0070085393 DJVM LLC
5908 SIOUX RIVER RD
SIOUX CITY, IA 51109
US
Shelf Life
45 days

Quality Plan Status
1. Do Not Use - Initial Draft

Comments

Related Quality Plans 1

☒ Test Methods ☐ Specifications ☐ Frequencies ☐ Locations ☐ Monitoring ☐ Sources ☐ Calculations ☐ Rules

Test Type All


Test Method	Certificate	Control Class	Active
-------------	-------------	---------------	--------



Chapter 6: Quality Plan - Edit

NEW History

7. To view the audit trail of changes to the quality plan information, open the desired quality plan and click on the 'History' button.




- a. The “Aggregate History” of changes will be displayed. Click on the ‘Item History’ hyperlink to view the item history.

Aggregate History is the audit trail of changes to the “Quality Plan Information” (the Active/Inactive checkbox, the Quality Plan Status, and the Comments) plus changes to the contents of the “Quality Plan Tabs.”

Item History is the audit trail of changes to only the “Quality Plan Information” (the Active/Inactive checkbox, the Quality Plan Status, and the Comments).

Note: Specifications are not included in this audit trail because they have their own audit trail application.



NEW Edit Quality Plan Information

8. To edit the quality plan information, open the desired quality plan, and click on the 'Edit' button.

QualityPatrol – Quality Plans & Specs

User Guide

Version 2025-10



Quality Plans > 703335 • MOR-P2-Cake5-Pkg • 0070085393

Quality Plan

703335 • MOR-P2-Cake5-Pkg • 0070085393

History Edit Clone Delete

- a. Edit the Active/Inactive checkbox, the Quality Plan Status or the Comments as desired, and click the 'Save Changes' button.

Note: The Product, Operation, and Customer are defined during quality plan created and cannot be changed thereafter.

Note: Shelf Life is defined at the product level.

Quality Plans > 703335 • MOR-P2-Cake5-Pkg • 0070085393

Quality Plan

703335 • MOR-P2-Cake5-Pkg • 0070085393

Active

Save Changes Cancel

Product
703335 • Cheesecake NY Strawberry

Operation
MOR-P2-Cake5-Pkg

Customer
0070085393 • DJVM LLC
5908 SIOUX RIVER RD
SIOUX CITY, IA 51109
US
Shelf Life
45 days

Quality Plan Status
1. Do Not Use - Initial Draft

Comments
This quality plan is for the NY Strawberry Cheesecake that is shipped to DJVM in Iowa.

NEW Clone Quality Plan

Cloning a quality plan is an efficient way to create a new quality plan based on an existing, similar quality plan. When cloning, the following data is copied exactly as configured in the quality plan being cloned:

- Quality plan comments.
- Attribute, VSI, and VMI test methods including the frequency, location, and CoA applicability.
- Alphanumeric inputs including the frequency, locations, and start-up criteria.
- Product and customer specifications.
- DashPatrol settings.
- Quality State rules and conditions.

9. To clone a quality plan, open the quality plan that most closely reflects the desired quality plan and click the 'Clone' button.

Quality Plans > 703335 • MOR-P2-Cake5-Pkg • 0070085393

Quality Plan

703335 • MOR-P2-Cake5-Pkg • 0070085393

Active

History Edit Clone Delete

QualityPatrol – Quality Plans & Specs

User Guide

Version 2025-10



- a. Click on the Product field.

Clone Quality Plan

Clone From (origin quality plan)	Clone To (destination quality plan)
Product 703335 • Cheesecake NY Strawberry	*Product None
Operation MOR-P2-Cake5-Pkg	*Operation MOR-P2-Cake5-Pkg
Customer 0070085393 • DJVM LLC 8321 WEST RD LOWELL, IN 46356 US	Customer None

Clone Cancel

- b. Use the Search and filters in the product picker to locate the desired product and then click on the desired product card.

Select Product

Enter product number or name...

Product Group Any Product Type Any

Status Any Type Any

1 2 3 4 5 6 7 ... 303 Page Go

12346 Active
Name Better Butter Crème Deluxe
Product Group TCCS BETTERCREME Test Group

- c. If the quality plan is for a specific customer, click on the Customer field, search for the desired customer in the customer picker, and then click on the desired customer card.

Clone Quality Plan

Clone From (origin quality plan)	Clone To (destination quality plan)
Product 703335 • Cheesecake NY Strawberry	*Product None
Operation MOR-P2-Cake5-Pkg	*Operation MOR-P2-Cake5-Pkg
Customer 0070085393 • DJVM LLC 8321 WEST RD LOWELL, IN 46356 US	Customer None

Clone Cancel

- d. Review the product and customer to confirm they are accurate and click on the 'Clone' button.

Clone Quality Plan

Clone From (origin quality plan)	Clone To (destination quality plan)
Product 703335 • Cheesecake NY Strawberry	*Product 12346 Butter Crème
Operation MOR-P2-Cake5-Pkg	*Operation MOR-P2-Cake5-Pkg
Customer 0070085393 • DJVM LLC 8321 WEST RD LOWELL, IN 46356 US	Customer 70007764 Esc. Sec. No.8 TORRES

Clone Cancel



- e. Click the 'Clone' button in the confirmation window.

Clone Quality Plan

To clone the origin quality plan to the destination quality plan, click the 'Clone' button. Otherwise, click the 'Cancel' button.

Product 12346 Better Butter Crème Deluxe

Operation MOR-P2-Cake5-Pkg

Clone

Cancel

Note: If restricted test methods are found that conflict with the new operation or workcenter, then a conflict message will appear advising that such test methods will be excluded from the clone.

NEW Delete Quality Plan

10. To delete a quality plan, open the desired quality plan, and click on the 'Delete' button and click on the 'Delete' button in the confirmation window.

Note: Quality plans cannot be deleted if the plan has test methods in it, or if a sample that uses the quality plan has been initiated. If either is the case and the quality plan is no longer in use, deactivate the plan and change the plan status as appropriate.

Quality Plans > 703335 MOR-P2-Cake5-Pkg 0070085393

Quality Plan

703335 MOR-P2-Cake5-Pkg 0070085393 Active

History

Edit

Clone

Delete



Chapter 7: Related Quality Plans

NEW View

11. To view related quality plans, click the “Related Quality Plans” text to expand the contents.

Note: The number of related quality plans is shown in the number icon.

Note: A related quality plan is a plan that has the same product, but a different operation or customer. Related plans are not user definable; they are determined by QualityPatrol logic.

Quality Plan

703335 • MOR-P2-Cake5-Pkg • 0070085393 Active

History Edit Clone Delete

Product: 703335 • Cheesecake NY Strawberry

Operation: MOR-P2-Cake5-Pkg

Customer: 0070085393 • DJVM LLC

5908 SIOUX RIVER RD

SIOUX CITY, IA 51109

US

Shelf Life: 45 days

Quality Plan Status: 1. Do Not Use - Initial Draft

Comments: This quality plan is for the NY Strawberry Cheesecake that is shipped to DJVM in Iowa.

Related Quality Plans 1

a. Click on the product number hyperlink to view a related quality plan.

Related Quality Plans 1

ARL-BIR

703335

MOR-P2-Cake5-Pkg


703335 This Plan




Chapter 8: Test Methods

NEW Add Test Methods


12. To add a test method to a quality plan, open the desired quality plan, click on the Test Methods tab, and click on the 'Add Test Method' button.



- a. Enter a test method number and click the 'Search' button, or use the desired filters to locate the desired test method.



- b. The "Results" tab will display the number of test methods that align with the search or filter criteria. Click on the "Results" tab, scroll to locate the desired test method, and click on the test method card.



- c. The selected test method will appear in the test methods data grid. To edit the test method settings, click on the 'Edit' button.

QualityPatrol – Quality Plans & Specs

User Guide

Version 2025-10



Test Method	Certificate	Control Class	Active
VS ARL-6200-Color-a Color (a)	<input type="checkbox"/>	Inherit [Standard]	<input checked="" type="checkbox"/>

d. Edit the Certificate, Control Class, and Active information as desired and click the 'Save Changes' button.

Note: To make the test method appear on certificates (CoAs), check the Certificate checkbox.

Note: Select the applicable Control Class.

- **Inherit [class]** means the test method is inheriting the class defined at the test property or test method. The actual class that is being inherited is shown in the [].
- **Not Considered** means the test method is not considered when determining the quality state. This is often called “run and record” because the test results, even if they are out of spec, will not result in a “Hold” or “Rejected” quality state.
- **Standard** means the results of the test method are considered when determining the quality state. If a test result for this test method is out of spec, the sample's quality state will be set to “Hold” or “Rejected” depending on the defined quality state rules.
- **Essential** means the test method is essential to the performance of the product and elevated scrutiny should be applied when determining the quality state of a sample with out of spec test results. If a test result for this test method is out of spec, the sample's quality state will be set to “Hold” or “Rejected” depending on the defined quality state rules.
- **Preventive** means that this test method is preventive control as defined by food safety regulations and is a core component of the facility's overall food safety plan. Elevated scrutiny should be applied when determining the quality state of a sample with out of specification test results. If a test result for this test method is out of spec, the sample's quality state will be set to “Hold” or “Rejected” depending on the defined quality state rules.

Note: Check the “Active” checkbox to make the test method active or inactive. Checked = Active. Unchecked = Inactive.

Test Method	Certificate	Control Class	Active
VS ARL-6200-Color-a Color (a)	<input type="checkbox"/>	Inherit [Standard]	<input checked="" type="checkbox"/>

13. To add another test method to the quality plan, click on the 'Add Test Method' button. Repeat the steps as outlined above.

Test Method	Certificate	Control Class	Active
VS ARL-6200-Color-a Color (a)	<input type="checkbox"/>	Inherit [Standard]	<input checked="" type="checkbox"/>

14. To remove a test method from a quality plan, click on the test method name dropdown and select

QualityPatrol – Quality Plans & Specs


User Guide

Version 2025-10




“Remove Test Method” from the submenu. Click the ‘Remove’ button in the confirmation window.

Note: Removing a test method from a plan has no effect on historical test results data. However, the subsequent QualityPatrol samples that are created will not include the test method as a required test.



15. To filter the test methods displayed, select the desired test type in the Test Type filter.










Chapter 9: Specifications

NEW Specification Levels

Specifications apply only to variable test methods. Alphanumeric and attribute data types do not have specifications. QualityPatrol supports specification management at the following levels:

- i.  **None:** No specifications have been defined.
- ii.  **Test Method Specifications:** For a given quality plan, the test method spec is used when: (i) a test method spec is defined; and (ii) no specification is defined at the product, customer, or ship-to address level.
- iii.  **Product Specifications:** For a given quality plan, the product spec is used when: (i) a product spec is defined; and (ii) no specification is defined at the customer or the ship-to address level. Accordingly, if both a test method spec and a product spec exist, the product spec is used.
- iv.  **Customer Specifications:** For a given quality plan, the customer spec is used when: (i) a customer spec is defined; and (ii) no specification is defined at the ship-to address level. Accordingly, if both a product spec and a customer spec exist, the customer spec is used.

At the customer level, specifications for the applicable customer can be further defined for the unique ship-to locations that are associated with the customer.

- v.  **Customer Ship-To Specifications:** For a given quality plan, the customer ship-to spec is used for all test plans that contain a test method with a defined customer ship-to spec, regardless of what other specifications exist.

NEW Specification Types

The following specification types can be defined:

- i. **Individual Specification:** An Individual Specification defines the acceptable tolerance for a single measurement taken from a process, product, or component.
 - Used when every unit must meet the specification on its own.
 - Ensures that no out-of-spec parts are produced or accepted.
- ii. **Average Specification:** An Average Specification sets the allowable range for the average value of multiple measurements — typically from a sample taken during a production run.
 - Common in statistical process control (SPC).
 - Allows for minor variation between individual parts as long as the process average is stable and within spec.



-
- iii. **Range Specification:** A Range Specification defines the maximum allowable difference between the largest and smallest values in a set of measurements.
- Focuses on consistency and process control.
 - High range values may indicate variation or instability in the process, even if averages are in spec.

NEW Specification Tolerance Limits

The following specification tolerance limits can be defined:

✂ **Upper Reasonable Data Limit (URDL):** The upper boundary of data values for a given process. Measurements beyond this limit are considered statistically unreasonable, potentially due to measurement error, data entry mistakes, or process failure.

☑ **Upper Specification Limit (USL):** The highest value a product characteristic or process parameter can have and still meet design or customer requirements.

⚠ **Upper Warning Limit (UWL):** A threshold set below the upper specification limit (USL) to provide an early warning that a process may be drifting out of control or approaching an unacceptable condition even though it's still within spec.

🎯 **Target:** The desired value that a measurement should achieve under optimal manufacturing conditions. It is typically the center point between the upper and lower specification limits and represents the goal for process or product performance.

⚠ **Lower Warning Limit (LWL):** A threshold set above the lower specification limit (LSL) to provide an early warning that a process may be drifting out of control or approaching an unacceptable condition even though it's still within spec.

☑ **Lower Specification Limit (LSL):** The lowest value a product characteristic or process parameter can have and still meet design or customer requirements.

✂ **Lower Reasonable Data Limit (LRDL):** The lower boundary of data values for a given process. Measurements beyond this limit are considered statistically unreasonable, potentially due to measurement error, data entry mistakes, or process failure.

NEW Test Method Specifications

🔧 **Test Method Specifications** Apply to test methods with a defined test method spec provided that a specification is not defined at the product, customer, or ship-to address level.

Use Case: If a test method always has the same specification regardless of the product or customer, then a test method level specification should be considered. This way, if the specification limits need to be edited, then this can be done in only one place (the test method spec) and the new spec will cascade to all test plans using the test method.

There are two ways to create a test method specification.


Method #1

16. From the PatrolSuite homepage, click on the Configure hyperlink in the QualityPatrol card and then click on the Test Methods card.


QualityPatrol – Quality Plans & Specs

User Guide


Version 2025-10




- a. Use the Search feature or the filters to locate the desired test method and click on the test method name hyperlink.



- b. Click on the “Specifications” tab then click on the ‘Specifications’ button.



- c. Click on the ‘Add New Spec’ button, or the ‘Edit’ button if the test method already has a specification.



- d. Enter the desired target, warning, specification, reasonable data limit values, and the Reason for Change, and click the ‘Save Changes’ button.

Note: To enter the upper and lower specs based on a calculation from target, enter the target, then enter either the +/- numerical value or the %, then click the ‘Calculate Individual Specs’ button.

QualityPatrol – Quality Plans & Specs

User Guide

Version 2025-10



Note: To enter the URDL and LRDL based on a calculation from target, enter the target then click on the 'Target x 2,' 'Target x 4,' 'Target x 9,' or the "100" button.

Method #2

17. From a quality plan, click on the desired test method dropdown and click on "View Test Method" submenu option.

a. Click on the "Specifications" tab then click on the 'Specifications' button.

b. Click on the 'Add New Spec' button, or the 'Edit' button if the test method already has a specification.



The screenshot shows two panels. The top panel is for 'Seconds' and has a red circle around the '+ Add New Spec' button. The bottom panel is for 'Feet per Minute' and has a red circle around the 'Edit' button. Below the 'Feet per Minute' header is a table with columns 'Individual', 'Average', and 'Range'. The table has two rows: 'Upper Reasonable Data Limit' and 'Upper Spec Limit'. The 'Upper Spec Limit' row shows a value of '60.0 Ft/Min'.

- c. Enter the desired target, warning, specification, reasonable data limit values, and the Reason for Change, and click the 'Save Changes' button.

Note: To enter the upper and lower specs based on a calculation from target, enter the target, then enter either the +/- numerical value or the %, then click the 'Calculate Individual Specs' button.

Note: To enter the URDL and LRDL based on a calculation from target, enter the target then click on the 'Target x 2,' 'Target x 4,' 'Target x 9,' or the "100" button.

The screenshot shows a 'Save Changes' dialog box. At the top right, there are 'Save Changes' and 'Cancel' buttons, with a red circle around the 'Save Changes' button. The dialog is divided into two main sections. The top section is titled 'Individual' and contains a table with columns for various limits and a 'Target' column. The table has rows for 'Upper Reasonable Data Limit', 'Upper Spec Limit', 'Upper Warning Limit', 'Target', 'Lower Warning Limit', 'Lower Spec Limit', and 'Lower Reasonable Data Limit'. The 'Target' row has a value of '100'. To the right of the table are buttons for 'Target x 2', 'Target x 4', 'Target x 9', and '100'. Below the table is a section for 'Reason for Change' with a text input field containing 'First specs issued by customer.' and a 'Calculate Individual Specs' button.

NEW Specifications Audit Trail

18. To view the audit trail of specification changes to a test method, open the desired test method, click on the 'Specifications' tab, and then click the 'Specifications' button.

QualityPatrol – Quality Plans & Specs

User Guide

Version 2025-10



Tests Methods > JDB-9701-Mix-Time

Test Method

VS JDB-9701-Mix-Time Active

Total mix time in seconds

General **Specifications** Monitoring Restrictions Required Files Test Costs

Test Method Specifications ⓘ

Seconds

History Clone Delete

a. Click on the “History” tab.

Test Method Specification ⓘ

Test Method JDB-9701-Mix-Time Total mix time in seconds

Product Any

Customer Any

Specifications History

b. The most recent revision will be displayed at the top. The spec value “From and To” changes are displayed in the Individual, Average, and Range columns as applicable.

Specifications History

Revision 2 Seconds Source QualityPatrol Chris Lloyd

	Individual	Average	Range
Upper Reasonable Data Limit	200		
Upper Spec Limit	From 110 To 120		
Upper Warning Limit	109		
Target	100		
Lower Warning Limit	91		
Lower Spec Limit	From 90 To 80		
Lower Reasonable Data Limit	20		

Reason for Change: Customer request #4412. Active Date 10/14/2025 4:39 PM | cor |

Revision 1 Seconds Source QualityPatrol Chris Lloyd

	Individual	Average	Range
Upper Reasonable Data Limit	200		



Chapter 10: Specification Editor

NEW Introduction

The Specification Editor enables the management of the specifications for test methods that are already part of a Quality Plan. Additional test methods cannot be added to a Quality Plan via the Specification Editor.


Specifications apply only to variable test methods. Alphanumeric and attribute data types do not have specifications. QualityPatrol

The specification level for each test method is identified with the following symbols.

- No specification
- Test method specification
- Product specification
- Product/Customer specification
- Product/Customer Ship-To specification

NEW View

19. To view all the specifications associated with a product quality plan, open the quality plan, and click on the 'Product Specifications' button.



- a. The specification editor will open, pre-filtered for the applicable product, operation, and customer. This page is organized into two sections: (i) Product information and filters; and (ii) The test methods which are organized by Test Property.



12345 > Specifications

PIZZA

12345 ▾ Nassau Pizza Thin Crust 12" Active

Product Type
FG

Test Group
None

Shelf Life
120 days

Product Description
Thin crust pizza exclusive to Nassau Pizza.

Data Type
Any

Test Type
Any

Customer
Any

Operation
Any

☒ Show only the test methods in the quality plan
☐ Show test methods from the Test Group quality plan

Reset

Temperature

Test Method	Specification	LRDL	LSL	LWL	Target	UWL	USL	URDL	Units
VS JDB-9800-Temperature-Mousse of Mousse	Temperature	Individual	15.2	38.0			43.0	68.8	Degrees F

Time

Test Method	Specification	LRDL	LSL	LWL	Target	UWL	USL	URDL	Units	
VS JDB-9701-Mix-Time	Total mix time in seconds	Individual	20	80	91	100	109	120	200	Seconds

Weights

Test Method	Specification	LRDL	LSL	LWL	Target	UWL	USL	URDL	Units
-------------	---------------	------	-----	-----	--------	-----	-----	------	-------

Note: Customer specifications are indented below the applicable test method. These specs are displayed as or because they technically apply to a product/customer combination.

VS MOR-0003-Volume-C	Height	Individual	40	50	60	N/A
Value at Point C						
0070007373 ▾ 344 Soria...		Individual	50	60	70	N/A
VS MOR-0003-Volume-D	Height	Individual	60	70	80	N/A
Value at Point D						
VS MOR-0003-Volume-Total		Individual	140	150	160	N/A
Calculated Sum of Height B+C+D						
70007373 ▾ 294 Soriana ...		Individual	145	155	165	N/A

NEW Define New Test Method Spec

Test Method Specifications: For a given quality plan, the test method spec is used when: (i) a test method spec is defined; and (ii) no specification is defined at the product, customer, or ship-to address level.

20. To define a new specification at the test method level, open a quality plan, click on the 'Product Specifications' button, and use the filters and the checkboxes to locate a test method of interest.

QualityPatrol – Quality Plans & Specs

User Guide

Version 2025-10



12345 > Specifications

PIZZA

12345 ▾ Nassau Pizza Thin Crust 12" Active Reset

Product Type
FG

Test Group
None

Shelf Life
120 days

Product Description
Thin crust pizza exclusive to Nassau Pizza.

Data Type Any ▾

Test Type Any ▾

Customer Any ▾

Operation Any ▾

☒ Show only the test methods in the quality plan
☐ Show test methods from the Test Group quality plan

- a. Click on the desired test method name with the icon and click on “Create Test Method Spec” in the submenu.

VM-3AA-3000-Diameter-X-Raw-Product-n3(5) ▾ Inches

Diameter
n3, 5 la Diameter X of Raw product by piece required samples n3, 5 lanes

Copy to Clipboard 60.00 80.00 100.00 Inches
60.00 80.00 100.00

VM View Test Method

Diameter
n3, 6 la Create Test Method Spec

- b. Click on the ‘Add New Spec’ button.

Specifications History

Seconds

Add New Spec

Individual

- c. Enter the desired target, warning, specification, and reasonable data limit values, record the Reason for Change, and click the ‘Save Changes’ button.

Note: To enter the upper and lower specs based on a calculation from target, enter the target, then enter either the +/- numerical value or the %, then click the ‘Calculate Individual Specs’ button.

Note: To enter the URDL and LRDL based on a calculation from target, enter the target then click on the ‘Target x 2,’ ‘Target x 4,’ ‘Target x 9,’ or the “100” button.

Specifications History

Seconds

Save Changes Cancel

Individual

Upper Reasonable Data Limit 200 Target x 2 Target x 4 Target x 9 100

Upper Spec Limit 110

Upper Warning Limit 109

Target 100 +/- 10

Lower Warning Limit 91 Calculate Individual Specs

Lower Spec Limit 90

Lower Reasonable Data Limit 20 Target x .5 Target x .25 Target x .1 0

Reason for Change
First specs issued by customer.



NEW Edit Existing Test Method Spec

Test Method Specifications: For a given quality plan, the test method spec is used when: (i) a test method spec is defined; and (ii) no specification is defined at the product, customer, or ship-to address level.

21. To edit an existing specification at the test method level, open a quality plan, click on the 'Product Specifications' button, and use the filters and the checkboxes to locate a test method of interest.

- a. Click on the desired test method name with the icon and click on "Edit Test Method Spec" in the submenu.

Test Method Name	Unit	Value 1	Value 2	Value 3	Value 4	Value 5	Value 6	Value 7
SAA-3000-Diameter-X-Raw-Product-n3(3)	Individual	10.00	100.00	120.00	150.00	180.00	200.00	500.00
Diameter n3, 3 la			100.00	120.00	150.00	180.00	200.00	
Diameter X of Raw product by piece required samples n3, 3 lanes			10.00				25.00	
Copy to Clipboard								Inches
View Test Method								
Edit Test Method Spec								
			60.00	80.00	100.00			Inches
			60.00	80.00	100.00			

- b. Click on the 'Edit' button.

- c. Edit the desired target, warning, specification, and reasonable data limit values, record the Reason for Change, and click the 'Save Changes' button.

Note: To enter the upper and lower specs based on a calculation from target, enter the target, then enter either the +/- numerical value or the %, then click the 'Calculate Individual Specs' button.

Note: To enter the URDL and LRDL based on a calculation from target, enter the target then click on the 'Target x 2,' 'Target x 4,' 'Target x 9,' or the "100' button.

Note: PatrolSuite can be integrated with an external system such as PLM software where specs are defined and managed. Specs defined and managed outside of QualityPatrol will not be editable from within QualityPatrol and will have an Externally Managed Data warning.



Specifications History

Seconds

Save Changes Cancel

Individual

Upper Reasonable Data Limit 200 Target x 2 Target x 4 Target x 9 100

Upper Spec Limit 110

Upper Warning Limit 109

Target 100 +/- 10

Lower Warning Limit 91

Lower Spec Limit 90

Lower Reasonable Data Limit 20 Target x .5 Target x .25 Target x .1 0

Calculate Individual Specs

Reason for Change

First specs issued by customer.

- d. To delete an existing test method specification, click the 'Delete' button and click the 'Delete' button in the confirmation window.

Specifications History

Inches

Edit Delete

	Individual	Average	Range
Upper Reasonable Data Limit	500.00 in		
Upper Spec Limit	200.00 in	200.00 in	25.00 in

NEW Define New Product Spec

Product Specifications: For a given quality plan, the product spec is used when: (i) a product spec is defined; and (ii) no specification is defined at the customer or the ship-to address level. Accordingly, if both a test method spec and a product spec exist, the product spec is used.

22. To define a new specification at the product level, open a quality plan, click on the 'Product Specifications' button, and use the filters and the checkboxes to locate a test method of interest.

12345 > Specifications

PIZZA

12345 ▾ Nassau Pizza Thin Crust 12" Active Reset

Product Type FG

Test Group None

Shelf Life 120 days

Product Description Thin crust pizza exclusive to Nassau Pizza.

Data Type Any

Test Type Any

Customer Any

Operation Any

☒ Show only the test methods in the quality plan

☐ Show test methods from the Test Group quality plan

- a. Click on the desired test method name and click on "Create Product Spec" in the submenu.

QualityPatrol – Quality Plans & Specs

User Guide

Version 2025-10



VS SAA-9801-PP-Humidity-Proof	% Humidity of	Percent Humidity
Proof		
VS % Humidity of Proofer	Number of belts	N/A
being		
VS Copy to Clipboard	Humidity of	Percent Humidity
Ambie		
VS View Test Method	rate of line	Minutes
VS Create Test Method Spec	Gap of	RPM
Sheet		
VS Create Product Spec		
VS Create Customer Spec		

b. Click on the 'Add New Spec' button.

Specifications History

Seconds

+ Add New Spec

Individual

c. Enter the desired target, warning, specification, and reasonable data limit values, record the Reason for Change, and click the 'Save Changes' button.

Note: To enter the upper and lower specs based on a calculation from target, enter the target, then enter either the +/- numerical value or the %, then click the 'Calculate Individual Specs' button.

Note: To enter the URDL and LRDL based on a calculation from target, enter the target then click on the 'Target x 2,' 'Target x 4,' 'Target x 9,' or the '100' button.

Specifications History

Percent Humidity

Save Changes Cancel

Individual

Upper Reasonable Data Limit

Upper Spec Limit 75.0

Upper Warning Limit

Target 70.0 +/- % Calculate Individual Specs

Lower Warning Limit

Lower Spec Limit 65.0

Lower Reasonable Data Limit

Target x .5 Target x .25 Target x .1 0

Reason for Change

Initial entry of customer specs.

d. The new specs will appear marked with the product icon.

Test Method	Specification	LRDL	LSL	LWL	Target	UWL	USL	URDL	Units
VS SAA-9801-PP-Humidity-Proof	% Humidity of	Individual	65.0		70.0		75.0		Percent Humidity
Proofer									

NEW Edit Existing Product Spec

Product Specifications: For a given quality plan, the product spec is used when: (i) a product spec is

QualityPatrol – Quality Plans & Specs


User Guide

Version 2025-10



defined; and (ii) no specification is defined at the customer or the ship-to address level. Accordingly, if both a test method spec and a product spec exist, the product spec is used.

23. To edit an existing specification at the product level, open a quality plan, click on the 'Product Specifications' button, and use the filters and the checkboxes to locate a test method of interest.

- a. Click on the desired test method name with the  icon and click on "Edit Product Spec" in the submenu.

- b. Click on the 'Edit' button.'

- c. Enter the desired target, warning, specification, and reasonable data limit values, record the Reason for Change, and click the 'Save Changes' button.

Note: To enter the upper and lower specs based on a calculation from target, enter the target, then enter either the +/- numerical value or the %, then click the 'Calculate Individual Specs' button.

Note: To enter the URDL and LRDL based on a calculation from target, enter the target then click on the 'Target x 2,' 'Target x 4,' 'Target x 9,' or the "100" button.

Note: PatrolSuite can be integrated with an external system such as PLM software where specs are defined and managed. Specs defined and managed outside of QualityPatrol will not be editable from within QualityPatrol and will have an Externally Managed Data warning.



Specifications History

Percent Humidity

Save Changes Cancel

Individual

Upper Reasonable Data Limit

Upper Spec Limit 75.0

Upper Warning Limit

Target 70.0 +/- %

Lower Warning Limit

Lower Spec Limit 65.0

Lower Reasonable Data Limit

Calculate Individual Specs

Reason for Change

Initial entry of customer specs.

NEW Define New Customer Spec

Customer Specifications: For a given quality plan, the customer spec is used when: (i) a customer spec is defined; and (ii) no specification is defined at the ship-to address level. Accordingly, if both a product spec and a customer spec exist, the customer spec is used.

24. To define a new specification at the customer level, open a quality plan, click on the 'Product Specifications' button, and use the filters and the checkboxes to locate a test method of interest.

12345 > Specifications

PIZZA

12345 Nassau Pizza Thin Crust 12" Active Reset

Product Type FG

Test Group None

Shelf Life 120 days

Product Description Thin crust pizza exclusive to Nassau Pizza.

Data Type Any

Test Type Any

Customer Any

Operation Any

☒ Show only the test methods in the quality plan

☐ Show test methods from the Test Group quality plan

- a. Click on the desired test method name and click on "Create Customer Spec" in the submenu.

VS SAA-9801-PP-Humidity-Proof	% Humidity of	Individual	65.0	70.0	75.0	Percent Humidity
Proof	% Humidity of Proofer					
VS being	Number of belts					N/A
VS Copy to Clipboard						
VS Ambie	Humidity of					Percent Humidity
VS View Test Method						
VS Edit Product Spec	ting rate of line					Minutes
VS Create Customer Spec	Gap of					RPM

- b. Select the desired customer from the business partner picker. Use the Search and filters as needed to locate the business partner.

QualityPatrol – Quality Plans & Specs

User Guide

Version 2025-10



Select Business Partner

70033921 Search

Country Any Market Any Priority Any Type Any

Cancel

70033921 Active Priority None Market

BON SUISSE INC FRF

14860 COMMUNITY RD STE 100

POWAY CA 92037 US

- c. Click on the 'Add New Spec' button.

Note: The applicable customer will be displayed.

- d. Enter the desired target, warning, specification, and reasonable data limit values, record the Reason for Change, and click the 'Save Changes' button.

Note: To enter the upper and lower specs based on a calculation from target, enter the target, then enter either the +/- numerical value or the %, then click the 'Calculate Individual Specs' button.

Note: To enter the URDL and LRDL based on a calculation from target, enter the target then click on the 'Target x 2,' 'Target x 4,' 'Target x 9,' or the '100' button.

QualityPatrol – Quality Plans & Specs

User Guide

Version 2025-10



The new specs will appear indented from the test method and marked with the applicable customer icon.

= Product/Customer

= Product/Customer Ship-To

VS SAA-9801-PP-Humidity-Proof	% Humidity of	Individual	65.0	70.0	75.0	Percent Humidity
Proofer						
70033921	BON SUISE INC FRF	Individual	70.0	75.0	80.0	Percent Humidity
70007373	294 Soriana Veracruz Floresta	Individual	145	155	165	N/A

NEW Edit Existing Customer Spec

Customer Specifications: For a given quality plan, the customer spec is used when: (i) a customer spec is defined; and (ii) no specification is defined at the ship-to address level. Accordingly, if both a product spec and a customer spec exist, the customer spec is used.

25. To edit an existing specification at the customer level, open a quality plan, click on the 'Product Specifications' button, and use the filters and the checkboxes to locate a test method of interest.

12345 > Specifications

PIZZA

12345 ▾ Nassau Pizza Thin Crust 12" Active Reset

Product Type
FG

Test Group
None

Shelf Life
120 days

Product Description
Thin crust pizza exclusive to Nassau Pizza.

Data Type Any ▾

Test Type Any ▾

Customer Any ▾

Operation Any ▾

☒ Show only the test methods in the quality plan

☐ Show test methods from the Test Group quality plan

- a. Click on the customer number dropdown for the desired test method name and click on "Edit Customer" in the submenu.

VS SAA-9801-PP-Humidity-Proof	% Humidity of	Individual	65.0	70.0	75.0	Percent Humidity
Proofer						
700339	BON SUI	Individual	70.0	75.0	80.0	Percent Humidity
11860	of belts					N/A
POWAY						
2	Individual	2	3	3		N/A
3	ity of					Percent Humidity
4						

- b. Click on the 'Edit' button.'

Customer Specification

Test Method SAA-9801-PP-Humidity-Proof % Humidity of Proofer

Product 12345 Nassau Pizza Thin Crust 12"

Customer 70033921 BON SUISE INC FRF

Specifications History

Percent Humidity Edit Delete

- c. Enter the desired target, warning, specification, and reasonable data limit values, record the Reason

QualityPatrol – Quality Plans & Specs

User Guide

Version 2025-10



for Change, and click the 'Save Changes' button.

Note: To enter the upper and lower specs based on a calculation from target, enter the target, then enter either the +/- numerical value or the %, then click the 'Calculate Individual Specs' button.

Note: To enter the URDL and LRDL based on a calculation from target, enter the target then click on the 'Target x 2,' 'Target x 4,' 'Target x 9,' or the "100" button.

Note: PatrolSuite can be integrated with an external system such as PLM software where specs are defined and managed. Specs defined and managed outside of QualityPatrol will not be editable from within QualityPatrol and will have an Externally Managed Data warning.

Specifications History

Percent Humidity

Save Changes Cancel

Individual

Upper Reasonable Data Limit Target x 2 Target x 4 Target x 9 100

Upper Spec Limit 75.0

Upper Warning Limit

Target 70.0 +/- %

Lower Warning Limit

Lower Spec Limit 65.0

Lower Reasonable Data Limit Target x .5 Target x .25 Target x .1 0

Reason for Change

Initial entry of customer specs.



Chapter 11: Search Specifications


NEW Introduction

The Search Specification page provides a broad view of specifications associated with any test method/product/customer combination including:


- i) Any test method regardless of the product or customer
- ii) Any product regardless of the test method or customer
- iii) Any customer regardless of the test method or product

NEW Search

26. To search for specifications, click the “Specifications” hyperlink in the QualityPatrol card.



- a. Use the Search feature or the filters to locate the desired specifications.





NEW Edit

27. To edit a specific specification for a test method/product/customer, click on the desired test method name and select “Edit Specification” from the submenu.

12

Test Method	Product	Customer	Specification	LRDL	LSL	LWL	Target	UWL	USL	URDL	Units
MUR-0002-Bake-Time	14138 Sealer ZTFCO 16oz		Individual				10.00				Minutes
Total Bake Time	14138 Sealer ZTFCO 16oz		Individual				325				Degrees F
Copy to Clipboard	14138 Sealer ZTFCO 16oz		Individual	0.00	11.10	11.50		11.90	21.42		Inches
View Test Method	14138 Sealer ZTFCO 16oz		Average		11.10	11.50		11.90			
	14138 Sealer ZTFCO 16oz		Range			0.00		0.80			
Edit Specification	14138 Sealer ZTFCO 16oz		Individual	0.00	11.10	11.50		11.90	21.42		Inches
	14138 Sealer ZTFCO 16oz		Average		11.10	11.50		11.90			

a. Alternatively, to open the Specification Editor for a desired product, click on the product number dropdown and select “Specifications” from the submenu.

MUR-0002-Bake-Time

14138

Individual

10.00

Minutes

MUR-0003-Oven-Temp-Convection

Product

14138

325

Degrees F

MUR-3000-Length-Fresh-L10-n3

Certificates

14138

0.00

11.10

11.50

11.90

21.42

Inches

MUR-3000-Length-Fresh-n3

Complaints

14138

11.10

11.50

11.90

Samples

14138

0.00

0.80

Quality Plans

14138

0.00

11.10

11.50

11.90

21.42

Inches

Specifications

14138

11.10

11.50

11.90



Chapter 12: Sample Frequencies


NEW Introduction

Within the context of a Quality Plan, sample frequencies define how often a sample is taken to measure product characteristics or process parameters. Sample frequencies are either “time based” (e.g. Every 30 minutes) or “event based” (i.e. Every time a product unit is produced).


The population of sample frequencies that can be consumed in any quality plan is defined on page: [{PatrolSuite URL}/Quality/Setup/SampleFrequencies](#)

NEW View

28. To view the frequencies associated with the test methods of a quality plan, open the quality plan and click on the “Frequencies” tab.



a. The existing frequencies and the Startup condition for each test method will be displayed.



NEW Startup

“Startup” is a sample type used in QualityPatrol to distinguish manufacturing “startup” samples from “production” samples and has implications for how such data is used in CoAs, DashPatrol, DataPatrol, and Metrics & Reports. **Startup is not a frequency. Do not create a frequency called “Startup.”** If this mistake has been made, then the Startup frequency should be deleted. If it cannot be deleted because it has been consumed in test plans or samples have been created, then it should be deactivated. Once deleted or deactivated, consider re-training the users regarding the difference between sample types (Production and Startup) and frequencies.



Nonetheless, whether or not a test method is performed on Startup samples is defined on the Frequency tab.

29. To select which test methods need to be completed on Startup samples, click the 'Edit' button.

Test Method	Frequency	Startup
AT ARL-6200-Color Color - Visual Standard	1 Hr	<input type="checkbox"/>
VS BVL-3000-Finished-Thickness Thickness of finished sample		<input type="checkbox"/>

- Check the checkbox(es) that correspond to the tests that need to be completed when Startup samples are taken and click the 'Save Changes' button.

Test Method	Frequency	Startup
AT ARL-6200-Color Color - Visual Standard	1 Hr	<input checked="" type="checkbox"/>
VS BVL-3000-Finished-Thickness Thickness of finished sample		<input checked="" type="checkbox"/>
VS BVL-3000-Finished-Width Width of finished sample		<input checked="" type="checkbox"/>
VS UNC-4901-Case-Weight Weight of 1 case		<input type="checkbox"/>

NEW Frequency Selection

30. To define how often samples need to be taken (and accordingly how often a test is performed), click on the 'Edit' button for the desired test method.

Test Method	Frequency	Startup
AT ARL-6200-Color Color - Visual Standard	1 Hr	<input checked="" type="checkbox"/>
VS BVL-3000-Finished-Thickness Thickness of finished sample		<input checked="" type="checkbox"/>

- In the Frequency Picker, use the scroll bar to locate the desired frequency, select the checkbox, and click the 'Save Changes' button.

Note: Checked = The frequency applies. Unchecked = The frequency does not apply.

Note: More than one frequency can be selected.



Select Frequencies

Automatic Frequencies

- ☒ 1 Hr
- ☐ 1/2 Batch
- ☐ 10 Min
- ☐ Each Unit
- ☐ Every other sample
- ☐ Prod Start

Manual Frequencies

- ☒ 1 Hr - Left Filler
- ☐ 1 Hr - Right Filler

- b. The updated quality plan will display the applicable frequencies and the applicability of the test method to Startup samples.

Test Methods Specifications Frequencies Locations Monitoring Sources Calculations Rules

Test Type All Product Specifications Add Test Method Edit

Test Method	Frequency	Startup
AT ARL-6200-Color Color - Visual Standard	1 Hr Prod End Re-Sample Speed Change	<input checked="" type="checkbox"/>

Edit



Chapter 13: Test Locations

NEW Introduction

Within the context of a Quality Plan, a test location is a physical location within the facility where tests are performed. For example, when a product sample is taken, 5 tests may need to be completed in the quality control lab, and 3 tests may need to be completed on the shop floor.


Additionally, since a product can be made at multiple facilities, and an operation can exist at multiple facilities, the quality plan for product XYZ may apply to multiple facilities. Accordingly, when the test location is selected, we need to identify the applicable facility. For example, Test Method 123 may be performed in the QC lab at facility #1 but is performed on the shop floor in Facility #2.

Test locations are defined for each facility on page: {PatrolSuite URL}/Manufacturing/TestLocations/{FacilityCode}

NEW View

31. To view the test locations for each test method within a quality plan, open the desired quality plan and click on the “Locations” tab.

Note: In the image below, Case Weight has no defined test location. In such cases, QualityPatrol will apply the facilities’ “Default” test location when a sample is initiated at such a facility.



NEW Edit

32. To edit the test locations for each test method within a quality plan, click on the ‘Edit’ button and select

QualityPatrol – Quality Plans & Specs

User Guide

Version 2025-10



the desired Facility in the submenu.

- Click the dropdown for the desired test method and select the test location where the test is performed. Then click the 'Save Changes' button.

- To define the test method test locations for another facility within a quality plan, click on the 'Edit' button and select the desired Facility in the submenu.

- Click the dropdown for the desired test method and select the test location where the test is performed. Then click the 'Save Changes' button.

- The applicable test locations will be displayed by facility.



Chapter 14: Monitoring

NEW Introduction

DashPatrol provides real-time visibility, via a Status Panel and live SPC charts, into the current manufacturing activities within a plant. DashPatrol is organized into two pages that provide content aimed at different audiences.


Facility Overview – The intended audience of the Facility Overview page is engineers and managers responsible for manufacturing and product quality who need to see the status of **only the most important test results**.

Line Details – The intended audience of the Line Details page is shop floor operators who need to see the status of **a broader set of Critical to Quality (CTQ) test results** for the purpose of maintaining process control and ensuring product quality.

The Monitoring tab defines which test methods within a test plan are displayed on the DashPatrol Facility Overview and Line Details pages.

NEW View

34. To view the monitoring settings for each test method within a quality plan, open the desired quality plan and click on the “Monitoring” tab.



Test Method	Facility Overview	Line Details
AT ARL-6200-Color Color - Visual Standard	<input checked="" type="checkbox"/>	<input type="checkbox"/>
VS BVL-3000-Finished-Thickness Thickness of finished sample	<input type="checkbox"/>	<input checked="" type="checkbox"/>

NEW Facility Overview

Note: The facility overview monitoring setting is controlled in the applicable test method. However, if the setting is changed here in the quality plan, then it is also changed in the test method, and that change



applies to every quality plan that contains that test method.

Note: As second criteria must be met for a test method to appear on the DashPatrol facility overview page:
The test method must also be checked as Line Details in the quality plan.

35. To edit the Facility Overview setting for a test method using the quality plan, open the desired quality plan and click on the “Monitoring” tab and then click on the ‘Edit’ button.

Test Method	Facility Overview	Line Details
AT ARL-6200-Color Color - Visual Standard	<input checked="" type="checkbox"/>	<input type="checkbox"/>
VS BVL-3000-Finished-Thickness Thickness of finished sample	<input type="checkbox"/>	<input checked="" type="checkbox"/>

- a. Check the desired checkbox and click the ‘Save Changes’ button. Remember that for a test method to appear in Facility Overview, the Line Details checkbox must also be checked.

Test Method	Facility Overview	Line Details
AT ARL-6200-Color Color - Visual Standard	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
VS BVL-3000-Finished-Thickness Thickness of finished sample	<input type="checkbox"/>	<input checked="" type="checkbox"/>

NEW Line Details

36. To edit the Line Details setting for a test method using the quality plan, open the desired quality plan and click on the “Monitoring” tab and then click on the ‘Edit’ button.

Test Method	Facility Overview	Line Details
AT ARL-6200-Color Color - Visual Standard	<input checked="" type="checkbox"/>	<input type="checkbox"/>
VS BVL-3000-Finished-Thickness Thickness of finished sample	<input type="checkbox"/>	<input checked="" type="checkbox"/>

- a. Check the desired Line Details checkbox and click the ‘Save Changes’ button.

Test Method	Facility Overview	Line Details
AT ARL-6200-Color Color - Visual Standard	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
VS BVL-3000-Finished-Thickness Thickness of finished sample	<input type="checkbox"/>	<input checked="" type="checkbox"/>



Chapter 15: Sources

NEW Introduction

The Sources tab identifies who requested that the test method be added to the quality plan. Additionally, whether the test appears on the certificate of analysis is also manageable on this tab.

Case Study: Quality plans inherently grow as people ask for additional testing for a product. Over the course of many years, test methods are added, but those added test are rarely removed. Accordingly, the Cost of Quality increases, often needlessly, as unnecessary testing is performed. It is a best practice to review and critique quality plans on a periodic basis. As part of this review, the test sources should be evaluated. As examples regarding the image below, one might question why the case weight is on the CoA if it is being measured at the request of engineering? Regarding the Vmag# requested by Engineering, does this test still need to be performed?

Test Method	Certificate	Test Required By
AT ARL-6200-Color Color - Visual Standard	<input checked="" type="checkbox"/>	Customer
VM BVL-2013-Vmag# Vmag number used at the Feeder	<input type="checkbox"/>	Engineering
VS BVL-3000-Finished-Thickness Thickness of finished sample	<input type="checkbox"/>	R&D
VS BVL-3000-Finished-Width Width of finished sample	<input type="checkbox"/>	R&D
VS UNC-4901-Case-Weight Weight of 1 case	<input checked="" type="checkbox"/>	Engineering

NEW View & Edit

37. To view the test sources for each test method within a quality plan, open the desired quality plan and click on the “Sources” tab.

Test Method	Certificate	Test Required By
AT ARL-6200-Color Color - Visual Standard	<input checked="" type="checkbox"/>	Customer
VM BVL-2013-Vmag# Vmag number used at the Feeder	<input type="checkbox"/>	Engineering
VS BVL-3000-Finished-Thickness Thickness of finished sample	<input type="checkbox"/>	R&D
VS BVL-3000-Finished-Width Width of finished sample	<input type="checkbox"/>	R&D
VS UNC-4901-Case-Weight Weight of 1 case	<input checked="" type="checkbox"/>	Engineering

- To edit the Source information, click the ‘Edit’ button, check the certificate checkbox to place the test method test results onto the CoA, select the Required By entity, and click the ‘Save Changes’ button.

QualityPatrol – Quality Plans & Specs

User Guide

Version 2025-10



Test Methods Specifications Frequencies Locations Monitoring Sources Calculations Rules

Test Type: All Product Specifications Save Changes Cancel

Test Method	Certificate	Test Required By
AT ARL-6200-Color Color - Visual Standard	<input checked="" type="checkbox"/>	Customer



Chapter 16: Calculations





NEW Introduction

Test methods can have “calculated results” based on other test method results **within the same test plan and within the same sample** (i.e. the frequencies have to match in order for the calculated test method to consume the “input” test methods). The quality plan identifies calculated test methods with the Sigma icon.

 VS DN-8801-Filling-Grade-C-%		Inherit [Standard]	<input checked="" type="checkbox"/>
Calculated C % Filling Grade			

NEW Create Expression

38. To create a calculated test method, first add all the test methods consumed in the calculation into the quality plan, then add the test method that reflects the calculation.

Test Methods Specifications Frequencies Locations Monitoring Sources Calculations Rules			
Test Type: All	Product Specifications	+ Add Test Method	Edit
Test Method	Certificate	Control Class	Active
 VS MOR-0003-Volume-B Height Value at Point B		Standard	<input checked="" type="checkbox"/>
 VS MOR-0003-Volume-C Height Value at Point C		Standard	<input checked="" type="checkbox"/>
 VS MOR-0003-Volume-D Height Value at Point D		Standard	<input checked="" type="checkbox"/>
 VS MOR-0003-Volume-Total Calculated Sum of Height B+C+D		Standard	<input checked="" type="checkbox"/>

- a. Click on the test method that reflects the calculation and select “Edit Expression” from the submenu. In this example, we are editing the test method that’s the result of calculating test methods Point B + Point C + Point D.

QualityPatrol – Quality Plans & Specs

User Guide

Version 2025-10



Test Methods			
Test Method	Certificate	Control Class	Active
VS MOR-0003-Volume-B Height Value at Point B		Standard	<input checked="" type="checkbox"/>
VS MOR-0003-Volume-C Height Value at Point C		Standard	<input checked="" type="checkbox"/>
VS MOR-0003-Volume-D Height Value at Point D		Standard	<input checked="" type="checkbox"/>
VS MOR-0003-Volume-Total Calculations		Standard	<input checked="" type="checkbox"/>

View Test Method
 Edit Expression
 Remove Test Method


b. The quality plan Expression Editor will open. Click on the Test Methods text to expand its contents.

Quality Plans > 703335 - MOR-P2-Cake5-Pkg - 0070085393 > Quality Plan > Edit Expression	
Edit Expression	
VS MOR-0003-Volume-Total	Test Methods Operators Functions Trigonometry Constants Other

c. Click on the first test method that comprises the calculation. The test method will be placed into the expression window on the left.

Quality Plans > 703335 - MOR-P2-Cake5-Pkg - 0070085393 > Quality Plan > Edit Expression	
Edit Expression	
VS MOR-0003-Volume-Total	Test Methods Test Method MOR-0003-Volume-B MOR-0003-Volume-C MOR-0003-Volume-D MOR-0003-Volume-Total Operators Functions Trigonometry Constants Other

d. Click on the “Operators” test to expand the mathematical operators. Click on the desired operator element on the right or simply click into the expression on the left and type the operator from your keyboard.

 Edit Expression

VS MOR-0003-Volume-Total

[MOR-0003-Volume-B] +

Save Changes

Cancel

▼ Test Methods

Test Method

MOR-0003-Volume-B

MOR-0003-Volume-C

MOR-0003-Volume-D


MOR-0003-Volume-Total

▼ Operators


Element	Description	Example
+	Add	1 + 1
-	Subtract	5 - 3
*	Multiply	2 * 2
/	Divide	12 / 3




- e. Choose the next test method (Volume C) add the “plus” operator element, and choose the last test method (Volume D) add click the ‘Save Changes’ button.



- f. On the Test Methods tab, the test method will be identified as “calculated.”




- g. On the Calculation stab, the calculation expression will be displayed.



NEW Test Expression

39. To test a calculation expression, click on the test method name and select “Edit Expression” from the submenu.



- a. Locate the Preview Expression area at the bottom of the Expression Editor. Click into each Test Value field and enter a numerical value and click on the ‘Calculate’ button. Evaluate the result to determine if the calculation produced the value that is expected. If not, click into the Edit Expression window, edit

QualityPatrol – Quality Plans & Specs

User Guide

Version 2025-10



the expression, and re-test.

40. There are numerous other highly technical expression possibilities, including if/then statements (e.g. “If the absolute value of the Auger Feed is greater than 2.0, then make the result of the TI-9400 attribute test method Pass. Otherwise, make the TI-9400 result Fail”). It would not be reasonable to document them all here. Contact MAP if you need assistance creating an expression

Note: In the example above, the if/then expression might look something like:

If(abs([AUGERFEED]) >2.0, 'Pass', 'Fail')



Chapter 17: Quality State Rules

NEW Introduction

The Rules tab enables the management of quality state rules that govern when and how a quality state is assigned to a sample for the applicable test plan. Additionally, these rules govern what is communicated to integrated business systems such as the ERP or MES. Note: For the communication to take place (i.e. QualityPatrol sending data to the local ERP or MES), the related integration work has to be performed.

A quality state rule consists of a:

Name. The name of the rule is user-defined and has no effect on logic.

Condition. The condition has to do with the test methods within a QualityPatrol sample, whether they are complete or not, and whether the result is within specifications or not. The conditions are not user definable. Contact MAP if a new condition will benefit your quality control process.

Assigned Quality State. This is the quality state that is assigned to the QualityPatrol sample if the conditions are met.

Notification. Once the quality state has been assigned, QualityPatrol can send a code to the financial system (e.g. ERP) or the manufacturing system (e.g. MES). The actual code that is sent is user definable in QualityPatrol/Configure/QualityStates. Contact MAP to request a quote to complete integration to enable this code to be sent to the ERP or MES.

QualityPatrol > Configure > Quality States

Quality States + New Quality State

Name	Type	Restricted	Approved with test results out of specification	Financial System Code	Manufacturing System Code
Approved: General Use	Approve			AGU	Approved

Note: One default quality state rule exists called “**All tests complete and all results in spec.**” This rule results in a quality state of “Approved” being assigned to a sample if:

- All tests are complete; and
- All results are within specifications.

Note: The above rule cannot be deleted and cannot be deactivated. However, the assigned quality state can be changed.



NEW View

41. To view the quality state rules associated a quality plan, open the quality plan and click on the “Rules” tab.

Condition	Assigned Quality State	Notify Financial System	Notify Manufacturing System	Active
All tests are complete and all results are within specification	Approved: General Use	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

NEW New Rule

42. To add a new rule, click on the ‘Add Rule’ button.

Condition	Assigned Quality State	Notify Financial System	Notify Manufacturing System	Active
All tests are complete and all results are within specification	Approved: General Use	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

a. Enter the rule Name, select the rule Condition and Assigned Quality State, check the Notifications checkboxes as desired, and click the ‘Save Changes’ button.

Quality Plans > 703335 - MOR-P2-Cake5-Pkg - 0070085393 > Quality State Rules > New Quality State Rule

New Quality State Rule Active

Name:

Condition: The results of any tests are out of specification

Assigned Quality State: Hold

☒ Notify Financial System

☒ Notify Manufacturing System

Note: Some of the Conditions are dependent on selected test methods in the quality plan. In such cases, select the desired test methods in addition to selecting the other setting noted above.

Condition: The results of any selected tests are out of specification

Assigned Quality State: Hold

☐ Notify Financial System

☒ Notify Manufacturing System

Test Methods

- ☒ VS MOR-0003-Volume-D Height Value at Point D
- ☒ VS MOR-0003-Volume-B Height Value at Point B
- ☒ VS MOR-0003-Volume-C Height Value at Point C
- ☐ VS MOR-0003-Volume-Total Calculated Sum of Height B+C+D

Note: Multiple quality state rules can be defined in one test plan.

QualityPatrol – Quality Plans & Specs

User Guide

Version 2025-10



Condition	Assigned Quality State	Notify Financial System	Notify Manufacturing System	Active
All tests are complete and all results are within specification	Approved: General Use	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
The results of any selected tests are out of specification	Hold	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

43. To edit a quality state rule, click on the desired rule name/condition.

Condition	Assigned Quality State	Notify Financial System	Notify Manufacturing System	Active
All tests are complete and all results are within specification	Approved: General Use	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
The results of any selected tests are out of specification	Hold	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

a. Click on the 'Edit' button.

Quality Plans > 703335 - MOR-P2-Cake5-Pkg - 0070085393 > Quality State Rules
 > The results of any selected tests are out of specification

Quality State Rule

The results of any selected tests are out of specification Active

Edit Delete

Assigned Quality State: Hold

☒ Notify Financial System
☒ Notify Manufacturing System

b. Check or uncheck the Active checkbox, make other setting changes as desired, and click on the 'Save Changes' button.

Quality Plans > 703335 - MOR-P2-Cake5-Pkg - 0070085393 > Quality State Rules > New Quality State Rule

New Quality State Rule Active

Save Changes Cancel

Name:

Condition: The results of any tests are out of specification

Assigned Quality State: Hold

☒ Notify Financial System
☒ Notify Manufacturing System

44. To delete a quality state rule, click on the desired rule name/condition.

Condition	Assigned Quality State	Notify Financial System	Notify Manufacturing System	Active
All tests are complete and all results are within specification	Approved: General Use	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
The results of any selected tests are out of specification	Hold	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

a. Click on the 'Delete' button and click on the 'Delete' button in the confirmation window.



Quality Plans > 703335 - MOR-P2-Cake5-Pkg - 0070085393 > Quality State Rules

> The results of any selected tests are out of specification

Quality State Rule

The results of any selected tests are out of specification

Active

Edit

Delete

Assigned Quality State

Hold

☒ Notify Financial System

☒ Notify Manufacturing System